



## Modular Liquid Handling Systems

HTS • Genomic Solutions • Protein Crystallography

**Gryphon**

**Optional 1-channel non-contact nano dispenser module:**  
- 100 nl up to 100 µl protein solution



**New! Optional LCP (Lipidic Cubic Phase) module:**  
- 25 nl to 2 µl of viscous solutions with a CV of < 5% in less than 90 seconds  
- Self-aligning syringe needle  
- Integrated LCP mixer  
- 96 well sandwich plate ready to seal in less than 2 minutes

**"All-In-One" Instruments:**  
Enable screen and protein drop setting in one protocol

**Phoenix**

8 plate positions and one wash station



**Optional up to 4-non-contact nano dispenser modules**

Also available from Art Robbins Instruments: [CrysCam-Microscope and IntelliPlates](#)

**Prima Xtallo**



- Preparation and reformatting of screening kits
- 2 pipettes (200 µl and 1000 µl)
- Platform for tubes and plates
- CV 0,2 - 2 % (precision of dispensing ± 1,5 % max.)



## The New D8 QUEST for Structural Biology

- 30% more intensity from the improved IµS™ with HELIOS MX optics
- New PHOTON 100 CMOS detector with three times larger active area
- High-performance, low-maintenance, all-air-cooled system, single phase power line
- Three year warranty on detector and X-ray source

Contact us for more details and system demonstration! [www.bruker.com/d8quest](http://www.bruker.com/d8quest)



Crystallography

Innovation with Integrity



Dunn Labortechnik GmbH · Thelenberg 6 · 53567 Asbach

Tel. +49 26 83 / 4 30 94 · Fax +49 26 83 / 4 27 76 · e-mail: [info@dunnlab.de](mailto:info@dunnlab.de) · Internet: [www.dunnlab.de](http://www.dunnlab.de)